CoCoALib - Feature #685

Use DivMask inside QBGenerator

16 Apr 2015 12:44 - John Abbott

Status:	New	Start date:	16 Apr 20)15	
Priority:	High	Due date:			
Assignee:		% Done:	0%		
Category:	Tidying	Estimated time:	0.00 hou)0 hour	
Target version:	CoCoALib-1.0	Spent time:	0.50 hou	r	
Description					
Internally QBGenera	tor usually makes very many div	isibility tests; use DivMask to make	these tests fa	ster.	
This will probably be	useful for current work on Implic	bitDrect			
Related issues:					
Related to CoCoALib - Bug #232: No test for QBGenerator			New	24 Sep 2012	
	Feature #775: PPMonoid with 1 extr		New	15 Sep 2015	

History

#1 - 16 Apr 2015 12:52 - John Abbott

We shall need collections of DivMask mirroring the elements of myCornerList and myAvoidList. A complication arises when myCornerList is sorted.

One possibility is to use PPWithMask inside myCornerList and myAvoidList, but this would force a change in the user interface (since both "lists" are publicly readable).

#2 - 16 Apr 2015 12:59 - John Abbott

- Description updated

#3 - 07 May 2015 14:36 - John Abbott

- Target version changed from CoCoALib-0.99536 June 2015 to CoCoALib-0.99540 Feb 2016

#4 - 02 Feb 2016 14:16 - John Abbott

- Target version changed from CoCoALib-0.99540 Feb 2016 to CoCoALib-1.0